

OtoClear® Ear Irrigation Tip

Bionix OtoClear® Ear Irrigation allows healthcare professionals to utilize a soft, flared, 3 streamered tip that eliminates the mess, improves patient safety, and is compatible with a variety of delivery devices to meet the personal preference of every facility.



PROBLEM

Traditional ear irrigation shoots water directly into the ear canal, which can be dangerous and run the risk of causing perforation to the eardrum.

SOLUTION

OtoClear® Ear Irrigation Tips direct water toward the ear canal wall, with three divergent streams, instead of directly at the Tympanic Membrane. The tip is also designed with a soft, flexible material for patient comfort and flared design to prevent backslash and over insertion.

RESULTS

A clean and safe procedure that improves efficiency and patient experience.

“Our office is very happy with the soft (OtoClear) tips. They are easier to maneuver in the ear canal and more comfortable for the patient.” - Laura Griffith, MA



Rx ONLY

OTOCLEAR TIPS ARE
MADE IN THE U.S.A.



Obelis s.a • Boulevard Général Wahis 53 • 1030 Brussels, BELGIUM
Tel: +(32) 2.732.59.54 • Fax: +(32) 2.732.60.03 • E-Mail: mail@obelis.net
Bionix Development Corp. • 5154 Enterprise Blvd. Toledo, OH, 43612, U.S.A.
Phone: +1 419.727.8421 • Fax: +1 800.455.5678 • www.Bionix.com



INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

800.551.7096 • www.bionix.com

©2019 Bionix Development Corp. | MKT-0040 | Rev. B

OtoClear® Ear Irrigation Tip

PRODUCT OFFERING



OtoClear SprayWash Kit

- Manual, intermittent lavage.
- Delivers greater volume compared to syringes or other manual devices.
- Efficient design allows for two-handed procedure.
- 500ml bottle volume.

#7290 - 1 SprayWash bottle, 20 OtoClear tips, 1 Ear Basin

#7295 - 1 SprayWash bottle, 5 OtoClear tips, 1 Ear Basin




OtoClear Aquabot® Kit

- Delivers greater volume compared to syringes or other manual devices.
- Efficient design allows for two-handed procedure.
- Provides continuous stream without manually pumping.
- 650ml bottle volume.

#7270 - 1 Aquabot® bottle, 20 OtoClear tips, 1 Ear Basin

#7275 - 1 Aquabot® bottle, 5 OtoClear tips, 1 Ear Basin



OtoClear Ear Irrigation Tips

- Directs water to the ear canal wall, with three divergent streams, instead of directly at the Tympanic Membrane.
- Flared design prevents backsplash and over insertion.

#7200 - OtoClear Ear Irrigation Tips (40/box)

#7202 - OtoClear Ear Irrigation Tips (20/box)

OtoClear Tip US Patent # 6,706,023.

Aquabot® is registered under Lunatech®, used under permission of Lunatech®, Inc. All rights reserved.



INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

800.551.7096 • www.bionix.com